

Trypsin / EDTA Solutions	Trypsin EDTA 10X in normal saline	Custom ATV
<b>Catalog Nos.</b>	<b>9342</b>	<b>9920</b>
<b>Component</b>	<b>mg/L</b>	<b>mg/L</b>
EDTA, Tetrasodium salt	—	226.17
Glucose	—	1000
Potassium Chloride	—	400
Sodium Chloride	8000	8000
Sodium Phosphate, dibasic Na <sub>2</sub> HPO <sub>4</sub> anhyd.	—	—
Phenol Red, Na salt	—	—
Trypsin 1:250	5000	500
EDTA, Disodium •2H <sub>2</sub> O	2214	—
Sodium Bicarbonate	—	580