

Trypsin / EDTA Solutions	Trypsin EDTA 10X in normal saline	Custom ATV
Catalog Nos.	9342	9920
Component	mg/L	mg/L
EDTA, Tetrasodium salt	—	226.17
Glucose	—	1000
Potassium Chloride	—	400
Sodium Chloride	8000	8000
Sodium Phosphate, dibasic Na ₂ HPO ₄ anhyd.	—	—
Phenol Red, Na salt	—	—
Trypsin 1:250	5000	500
EDTA, Disodium •2H ₂ O	2214	—
Sodium Bicarbonate	—	580