

Trypsin EDTA Solutions

Trypsin/EDTA Solutions	Versene solution	Trypsin 10x in normal solution	Trypsin EDTA 1x in normal saline	Trypsin EDTA 1x in HBSS
Catalog No.	9314	9336	9340	9341
Component	mg/L	mg/L	mg/L	mg/L
EDTA, Tetrasodium salt	200		226.17	226.17
Glucose (Dextrose) anhyd				1000
Potassium Chloride	200			400
Potassium Phosphate monobasic KH_2PO_4	200			60
Sodium Chloride	8000	5875	8768	8000
Sodium Phosphate, dibasic Na_2HPO_4 anhyd	1150			47.68
Phenol red, Na salt			10	10
Trypsin, 1:250		25000	500	500