

# Heavy Oil for Embryo Culture

An extra blanket of protection

Heavy Oil for Embryo Culture is a sterile, heavy, mineral oil used as an overlay when culturing small volumes of media to protect the fertilized embryo - preventing evaporation, changes in osmolality, and pH shifts, while providing the optimal *in vitro* environment for development.

Providing an extra blanket of protection for your embryo culture with optimal handling viscosity, each lot of Heavy Oil is validated by rigorous quality testing procedures for sterility, endotoxin, HSSA and MEA, with POV values reported on the COA for clear visibility.

- Ideal handling viscosity for IVF processes
- Superior quality control with MEGA - Our unique genetic mouse embryo assay for added security
- POV values reported on the COA to assure quality and confidence

We  
Have You  
Covered



For more information contact your local sales representative or visit [www.irvinesci.com](http://www.irvinesci.com).

## Ordering Information

Media	Catalog #	Size*	Additional Information	Shelf Life*	Storage
Heavy Oil for Embryo Culture	90189	100 mL 500 mL	Ready-to-use sterile heavy mineral oil overlay for small media volumes	2 years or 8 weeks after opening	Storage at 2-8°C

\*From date of manufacture