

Transfectory Supplement

Catalog #	Product	Size
91148	Transfectory Supplement	500mL

Intended Use

For research or further manufacturing use only.

Product Description

Transfectory Supplement is concentrated hydrolysate supplement designed to enhance cell growth and protein expression of CHO and various cell types. This supplement is free of animal-derived components.

Quality Assurance

All quality control test results are reported on a lot specific Certificate of Analysis which is available upon request.

Storage Instructions and Stability

Handle using aseptic techniques to avoid contamination. Store at 2-8°C away from light. Do not use after the assigned expiration date. Do not use any bottle of medium that shows evidence of particulate matter or cloudiness.

Precautions

The safety and efficacy of this product in diagnostic or other clinical uses has not been established. Please refer to the Safety Data Sheet for information regarding hazards and safe handling practices.

Directions for Use

For Stable Cell Lines

Supplement 10% of the initial culture volume (e.g. 3mL for 30mL culture in a shake flask) when cell density reaches approximately one third of the maximum viable cell density.

For Post-Transient-Transfection Use

Supplement 10% of the initial culture volume (e.g. 3mL for 30mL culture in a shake flask) one day after transfection to boost cell growth and titer.

Note: This supplement may interfere with transfection that uses cation transfection reagents.

Technical Support

CONTACT US

For more information or assistance contact Customer Service at:

- Email: fisitmrequest@fujifilm.com
- Direct line: +1 800 577 6097

WEBSITE RESOURCES

Visit the website at www.irvinesci.com for technical resources and information including:

- Safety Data Sheets (SDS)
- Certificate of Analysis (CoA) (when available)
- FAQs
- Product literature
- Complete list of offices and contact information by country

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705-5588 USA

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • www.irvinesci.com

© 2020 FUJIFILM Irvine Scientific, Inc. All rights reserved. FUJIFILM Irvine Scientific and its logo are registered trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. P/N 41061, Rev.01