


MAKLER[®] INSEMINATION DEVICE



A device to keep the Makler Cannula in place during the insemination, preventing leakage of the injected sperm.



SEFI-MEDICAL INSTRUMENTS LTD.

Distributed By: Irvine Scientific 2511 Daimler St. Santa Ana, CA 92705
TEL: 800-437-5706 FAX: 949-261-6522 www.irvinesci.com

Distributed By



IrvineScientific[®]

Grow With Us

MAKLER[®] INSEMINATION DEVICE

UNIQUE FEATURES AND ADVANTAGES

The Makler Insemination Device serves as a syringe holder to keep the Makler Cannula in place, preventing leakage of the injected sperm.

EASY MANIPULATION -

Due to its length, it is easy to introduce the Makler Cannula into the uterine cavity.



SEALING MECHANISM -

Prevents leakage of the injected sperm. Provides excellent sealing by pressing the Makler Cannula against the cervix.

ADJUSTABLE -

The pressure on the cervix can be adjusted by a clamp attached to the speculum.

KEEPS THE CANNULA IN PLACE- No need to hold it during the insemination.



SEFI MEDICAL INSTRUMENTS LTD.

Distributed By: Irvine Scientific 2511 Daimler St. Santa Ana, CA 92705
TEL: 800-437-5706 FAX: 949-261-6522 www.irvinesci.com

Distributed By



IrvineScientific
Grow With Us