

SAFETY DATA SHEET

Issuing Date: 24-Oct-2022

Version 1

CTGrade rh IL-2 C126S

1. PRODUCT AND COMPANY IDENTIFICATION

Product Name CTGrade rh IL-2 C126S

Product code 500-01

Product Use Laboratory chemicals.

Manufactured by
 FUJIFILM Irvine Scientific, Inc.
 2511 Daimler Street
 Santa Ana, CA 92705
 U.S.A.

SDSs are available at the following website(s): <http://www.irvinesci.com>

Company Phone Number Technical: 800-437-5706
 Fax: 949-261-6522

Emergency Telephone Transport-CHEMTREC Inside NA: 800-424-9300
 Transport-CHEMTREC Outside NA: 703-527-3887
 Transport-CANUTEC Inside Canada: 613-996-6666
 Medical Emergency (24 hour): 877-935-7387

E-mail fisitmrequest@fujifilm.com

2. HAZARDS IDENTIFICATION

Classification

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS).

GHS Label elements, including precautionary statements

Not classified

Precautionary Statements**Prevention**
Not applicable**Response**
Not applicable**Storage**
Not applicable**Disposal**
Not applicable**Hazards not otherwise classified (HNOC)**
Not classified**3. COMPOSITION/INFORMATION ON INGREDIENTS**

The product contains no substances which at their given concentration, are considered to be hazardous to health.

4. FIRST AID MEASURES**Description of first aid measures**

General advice	If symptoms persist, call a physician.
Eye contact	Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.
Skin contact	Wash off immediately with plenty of water. Get medical attention if irritation develops and persists.
Inhalation	Move to fresh air. If symptoms persist, call a physician.
Ingestion	If swallowed, do not induce vomiting - seek medical advice.

Most important symptoms and effects, both acute and delayed

None known.

Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

5. FIRE-FIGHTING MEASURES**Suitable Extinguishing Media**

Dry chemical, CO2, water spray or alcohol-resistant foam.

Unsuitable Extinguishing Media

None known.

Specific hazards arising from the chemical

None known.

Explosion Data

Sensitivity to Mechanical Impact none

Sensitivity to Static Discharge none

Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES**Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment. Avoid contact with the skin and the eyes.

Environmental precautions

Prevent product from entering drains.

Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Use personal protective equipment. Soak up with inert absorbent material. Pick up and transfer to properly labeled containers. Clean contaminated surface thoroughly. After cleaning, flush away traces with water.

7. HANDLING AND STORAGE**Precautions for safe handling**

Wear personal protective equipment. Avoid contact with skin, eyes or clothing. Avoid breathing vapors or mists. Ensure adequate ventilation.

Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place. Keep at temperatures below -20 °C / -4 °F.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION**Control parameters**

Exposure Guidelines This product does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Appropriate engineering controls

Engineering Measures Ventilation systems

Individual protection measures, such as personal protective equipment

Eye/Face Protection If splashes are likely to occur, wear safety glasses with side-shields.

Skin and body protection Wear protective gloves/clothing.

Respiratory protection	None required under normal usage. If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn.
General Hygiene Considerations	When using do not eat, drink or smoke. Take off contaminated clothing and wash before reuse. Regular cleaning of equipment, work area and clothing is recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	White	Physical state	Solid
Form	lyophilisate, Powder		
Odor	Not available	Odor Threshold	Not available
pH	No data available		
Specific Gravity	Not available	Molecular Weight	Not available
Flash point	Not available	Autoignition temperature	Not available
Decomposition temperature	Not available	Boiling point / boiling range	Not available
Melting point / melting range	Not available	Freezing Point	Not available
Flammability Limit in Air	Not available		
Oxidizing Properties	Not available	Explosive Properties	Not available
Solubility	Not available	Partition coefficient	Not available
Evaporation rate	Not available	Vapor Pressure	Not available
Vapor density	No data available	Density	Not available
VOC (lb/gal)	Not available	VOC (g/l)	Not available
Dynamic viscosity	Not available		

10. STABILITY AND REACTIVITY

Reactivity

Stable under recommended storage conditions.

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

None under normal processing.

Conditions to Avoid

Protect from light. Temperatures above -20 °C / -4 °F.

Incompatible materials

Strong oxidizing agents. Strong acids. Strong bases.

Hazardous decomposition products

None known.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information**Acute toxicity****Inhalation****Eyes****Skin****Ingestion**

Non-irritating during normal use.

Non-irritating during normal use.

Non-irritating during normal use.

Ingestion may cause stomach discomfort.

Component Information

No information available

Symptoms related to the physical, chemical and toxicological characteristics

No information available.

Delayed and immediate effects as well as chronic effects from short and long-term exposure**Irritation**

Non-irritating during normal use.

Corrosivity

No information available.

Sensitization

No information available.

Mutagenic Effects

No information available.

Reproductive Toxicity

No information available.

Carcinogenicity

None known.

STOT - single exposure

No information available.

STOT - repeated exposure

No information available.

Target organ effects

None known.

Aspiration hazard

No information available.

Numerical measures of toxicity - Product Information**ATEmix (oral)**

>5000 mg/kg

ATE: Acute toxicity estimate**12. ECOLOGICAL INFORMATION****Ecotoxicity**

The environmental impact of this product has not been fully investigated.

Persistence and degradability

No information available.

Bioaccumulation

No information available

Mobility

No information available.

Other adverse effects

No information available.

13. DISPOSAL CONSIDERATIONS

Waste Disposal Methods

Dispose of in accordance with local regulations.

Contaminated packaging

Do not re-use empty containers.

14. TRANSPORT INFORMATION

<u>DOT</u>	Not regulated
<u>TDG</u>	Not regulated
<u>MEX</u>	Not regulated
<u>ICAO (air)</u>	Not regulated
<u>IATA</u>	Not regulated
<u>IMDG</u>	Not regulated
<u>ADR</u>	Not regulated
<u>ADN</u>	Not regulated

15. REGULATORY INFORMATION

International Inventories

TSCA	No
DSL/NDSL	No
PICCS	No
EINECS/ELINCS	No
ENCS	No
IECSC	No
KECL	No
AIICS	No

***Yes - All component(s) of this product are included or are exempt from listing on the inventory.**

***No - Indicates the component(s) of this product are either not listed or have not been determined to be listed on the inventory.**

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AIICS - Australian Inventory of Industrial Chemicals

U.S. Federal Regulations

TSCA Sections 4, 5 and 12(b)

This product does not contain any chemicals regulated by TSCA Sections 4, 5 or 12(b).

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

SARA 311/312 Hazard Categories

Classification is shown in section 2 of this SDS.

Clean Water Act.

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

U.S. State Regulations**California Proposition 65**

This product does not contain any Proposition 65 chemicals

U.S. State Right-to-Know Regulations

This product does not contain any substances regulated by state right-to-know regulations.

International Regulations**Canada - NDSL**

This product does not contain any NDSL chemicals.

Mexico - Grade

No information available

Mexico - Carcinogen Status and Exposure Limits

No information available

Other Regulations

No information available.

16. OTHER INFORMATION

NFPA	Health hazards 1	Flammability 0	Instability 0	Special hazards -
HMIS	Health hazards 1	Flammability 0	Physical hazards 0	Personal protection B

Prepared By FUJIFILM Environment, Health and Safety, phone: 800-473-3854

Revision Date 24-Oct-2022

Revision Note No information available

Disclaimer The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

end